

# AUTOMATED SURFACE SWAB SAMPLING: A STATISTICAL COMPARISON OF A NOVEL APPROACH TO EXISTING METHODS.


PRESENTED BY:

KEITH BADER,

VP CLEANING SCIENCE & LABORATORY SERVICES

HYDE ENGINEERING + CONSULTING, INC.



28 April 2025



# Direct Sampling



- Regulatory Health Authorities suggest the use of direct and indirect samples for cleaning validation
  - Health Canada, Cleaning Validation Guide [1]
  - World Health Organization, Good Manufacturing Practices, Appendix 3, Cleaning Validation [2]
  - US FDA, Validation of Cleaning Processes [3]
  - Eudralex Volume 4, Annex 14 [4]
- While it's possible, depending on product and process characteristics, to rely on rinse sampling, most regulators expect surface swab sampling



# Sampling Considerations




Confined Space and Manual Swabbing  
vs. Remote Sampling Methods

- Sampling Sites Often Difficult to Reach and Require Confined Space Entry
- Confined Space Entry Increases Safety Risks as well as the Possibility of Equipment Damage
- Remote Devices Can be Validated and Used for Sample Collection



# Confined Space Entry Considerations

- Fall Protection
- Crane
- Entry Ladders
- Atmospheric Monitoring
- Forced Ventilation
- Additional Attendants
- EMS Personnel
- Confined Space Entry Increases the Likelihood of Sample Contamination





# Industry Approach to Remote Sampling



Recovery should be shown to be possible from all product contact materials sampled in the equipment with all the sampling methods used.



# Development and Qualification of Remote Method




- Control of the Swab Extension Pole is Difficult
- Inconsistent Coverage and Sampling Pattern
- Variable Pressure on the Swab Head can Impact Recovery
- Training, Practice, and Requalification are of Paramount Importance

# Study Design



- Hyde compared three different swab methods
  - Manual or hand swabbing
  - Swabbing with an extension pole
  - Swabbot's prototype automated swabbing device
- To ensure that the solutions used for the recovery performance characterization study are reliable
  - Carbon content characterization was performed for sucrose and bovine serum albumin solutions.
  - Three concentrations as well as blanks samples were analyzed.



# Study – Materials and Methods



## Sucrose (ACS Grade)

Blank (4x)  
0.5  $\mu\text{g}/\text{cm}^2$  (4x)  
1  $\mu\text{g}/\text{cm}^2$  (4x)  
5  $\mu\text{g}/\text{cm}^2$  (4x)



## BSA (1 mg/mL)

Blank (4x)  
0.5  $\mu\text{g}/\text{cm}^2$  (4x)  
1  $\mu\text{g}/\text{cm}^2$  (4x)  
5  $\mu\text{g}/\text{cm}^2$  (4x)



## Coupons

316L Stainless Steel  
20 Ra Surface Finish  
2.5" x 2.5"



## Total Organic Carbon

Sievers M9  
Acid Flowrate (1.0  $\mu\text{L}/\text{min}$ )  
Oxidizer Flowrate (1.0  $\mu\text{L}/\text{min}$ )




## Data Analysis

% Recovery  
• Blank adjusted /  
Based on positive  
controls  
One-Way ANOVA






# Swab Patterns



## Manual and Automated Swabbing



## Remote Swabbing



# Remote Swab Sampling



# Study Results




## Swab Sampling Results Summary (316L Stainless Steel)


Expected Concentration (ppm C)	Method	Average Recovery (%)		Average Recovery % SD		% Difference from Swabbot	
		Sucrose	BSA	Sucrose	BSA	Sucrose	BSA
0.5	Hand	100	107	1	3	1	6
1	Hand	95	102	7	6	3	3
5	Hand	96	101	7	2	3	5
0.5	Remote	89	90	3	4	11	11
1	Remote	90	86	5	10	8	14
5	Remote	92	82	5	8	7	16
0.5	Swabbot	99	101	1	3		
1	Swabbot	98	99	1	4		
5	Swabbot	99	96	1	6		

# Study Results

Sucrose Means Comparison Chart  
Red intervals that do not overlap differ.



BSA Means Comparison Chart  
Red intervals that do not overlap differ.








# Using Recovery Data to Improve Instrument Design

A comparison of swabbing patterns for data-driven instrument footprint reduction

# Swab Pattern Comparison



- Initial Prototype Design used the Pattern Shown Below
- Rotation of the Swab to a Perpendicular Orientation Required More Components and Operational Space
- Determine if the Pattern Used for Swab Extension Poles is Comparable to Rotation




# Section Break



**316L SS Swabbot Sampling Method Summary Table**

<b>Expected Concentration (ppm C)</b>	<b>Sampling Method</b>	<b>Avg Corr TOC (ppm C)</b>	<b>Recovery</b>	<b>Repeatability Recovery SD</b>
0.5	Original	0.453	97%	3%
1	Original	0.907	93%	2%
5	Original	4.77	91%	2%
0	Stable Orientation	0.0824	N/A	N/A
0.5	Stable Orientation	0.477	102%	1%
1	Stable Orientation	0.928	95%	1%
5	Stable Orientation	4.86	93%	1%
0	No Flip	0.0834	N/A	N/A
0.5	No Flip	0.456	98%	3%
1	No Flip	0.850	87%	1%
5	No Flip	4.63	89%	3%


# Resulting Design



- Cleanable enclosure
- Silicone Bumpers as Point of Contact with Equipment to Prevent Damage
- Resulting Design Weighs Less than 3 Kilograms
- The instrument is mounted on a carbon fiber extension pole.
- The instrument is deployed to the sample site using an equipment access system designed to facilitate positioning and prevent contamination or equipment damage.



# Deployment Examples



# Acknowledgements



Rick Mineo

- Founder, Swabbot, Inc.



Nicole Collier

- Principal Engineer, Hyde Engineering



Michael Lund

- Senior Scientist, Hyde Engineering



# References



1. Health Canada's Cleaning validation guide (GUI-0028), <https://www.canada.ca/en/health-canada/services/drugs-health-products/compliance-enforcement/good-manufacturing-practices/validation/cleaning-validation-guidelines-guide-0028/document.html#s9-1>, accessed November 14, 2022
2. The FDA Guide to Inspections : Validation of Cleaning Processes states <https://www.fda.gov/inspections-compliance-enforcement-and-criminal-investigations/inspection-guides/validation-cleaning-processes-793>
3. WHO TRS 1019 - Annex 3: Good manufacturing practices: guidelines on validation, Appendix 3 [https://cdn.who.int/media/docs/default-source/medicines/norms-and-standards/guidelines/production/who-good-manufacturing-practices-guidelines-on-validation.pdf?sfvrsn=9440a5c\\_0&download=true](https://cdn.who.int/media/docs/default-source/medicines/norms-and-standards/guidelines/production/who-good-manufacturing-practices-guidelines-on-validation.pdf?sfvrsn=9440a5c_0&download=true), accessed November 14, 2022
4. PICS Annex 15 " in section 10.12, <https://picscheme.org/docview/4590>